

MELODY



MELODY es una silla de baño en aluminio, para niños con dificultades motoras. Permite posicionar al niño con más comodidad y seguridad. Articulación del respaldo y piernas. Respaldo en red con varios colores disponibles (azul, blanco, verde, coral, rosa, ...). Materiales resistentes a limpieza y desinfección. Varios accesorios disponibles (apoyos laterales de cabeza, apoyo de tronco, cinturón de piernas para control de la abducción y aducción, ...).

MELODY is a bath chair with aluminium structure for kids with motor disabilities. It allows the positioning with comfort and security. Backrest and legs articulation. Comfortable mesh back easy to clean, available in several colors (blue, white, green, coral, pink, ...). Cleaning and disinfection resistant materials. Several accessories available (lateral head supports, trunk support, leg belts controlling the abduction and adduction, ...).

ORTHOSXXI, Unipessoal Lda.

Rua Santa Leocádia nº2735 // 4809-012 Sta. Leocádia de Briteiros // Guimarães // Portugal

Phone + 351 253 470 276 // Mobile +351 936 199 005 // Email margarida.silva@orthosxxi.com // Website www.orthosxxi.com

Características / Characteristics



Con base antideslizante para mayor seguridad

Anti-slip base for added safety



Red confortable disponible en varios colores

Comfortable mesh available in several colors



Regulación del respaldo y piernas

Backrest and legrest recline



MELODY			
Marca/ Brand: Orthos XXI			
Ancho asiento <i>Seat width</i>	310mm	Ancho total <i>Total width</i>	397mm
Profundidad asiento <i>Seat depth</i>	220mm	Longitud total <i>Total length</i>	550mm
Altura respaldo <i>Back height</i>	525mm	Altura máxima usuario <i>Maximum user height</i>	1000mm
Inclinación asiento <i>Seat incline</i>	4°	Peso máximo usuario <i>Maximum user weight</i>	25Kg

Garantía/ Warranty

2 años contra defectos de fabricación
2 years against manufacturing defects

Tiempo estimado de entrega / Delivery time

Consultar/ Consult

ORTHOSXXI, Unipessoal Lda.

Rua Santa Leocádia nº2735 // 4809-012 Sta. Leocádia de Briteiros // Guimarães // Portugal

Phone + 351 253 470 276 // Mobile +351 936 199 005 // Email margarida.silva@orthosxxi.com // Website www.orthosxxi.com